BLS	OT	AM	FR	EMR	AFA	MAFA
	✓		✓	✓	✓	✓

Steps

Ensure you are always wearing proper PPE before performing the skill.

Set-Up:

- 1. Ensure oxygen cylinder is safely handled and positioned.
- 2. Check the cylinder to confirm that it contains medical-grade (USP) oxygen.
- 3. Clear valve:
 - a) Ensure that you are at a safe distance from others and that the valve is facing away from yourself and others caution anyone nearby of your intent.
 - b) Remove protective covering and set aside O-ring.
 - c) Quickly open the cylinder's valve (for a maximum of one second) to clear debris.
- 4. Attach pressure regulator:
 - a) Confirm that the gasket (O-ring) is properly inserted/positioned.
 - b) Check to see that the pin index corresponds to the oxygen tank.
 - c) Seat the two metal prongs of the regulator inside the holes of the valve stem.
 - d) Hand-tighten the screw until the regulator is snug.
 - e) Open the cylinder by giving the valve one full turn and listen for leaks.
 - f) Examine the regulator's gauge to determine the pressure in the cylinder.

ADMINISTERING OXYGEN

- 1. Attach the delivery device to the oxygen port on the regulator.
- 2. Set the appropriate flow rate.
- 3. Listen and feel to make sure that oxygen is flowing into the delivery device.
- 4. If patient is responsive, explain the procedure.
- 5. Place the delivery device on the patient.

Organizational/local protocol for:						

Skill Notes

- If using a delivery device with an oxygen reservoir bag, ensure that it is full.
- Regulators and cylinders are pin indexed, meaning a regulator will fit into only the cylinder it was designed for. Be sure to confirm
 that the pin index on the regulator corresponds to that on the oxygen cylinder.
- To safely administer oxygen from an oxygen cylinder, you must first reduce the pressure of the oxygen inside from 2,000 pounds per square in. (psi) to a range of 30 to 70 psi. This is done using the oxygen regulator, which is either integrated or independent.